

Apex Resultsviewer

Javascript notation

Most programming languages

```
int myVariable = 5;
```

```
string myVariable = "blabla";
```

```
ObjectType myVariable = new ObjectType();
```

Javascript

```
var myVariable = whateveryouwant;
```

```
myVariable = Object;  
myVariable.myNumber = 5;  
myVariable["myNumber"] = 5;  
myVariable.myString = "bla";  
myVariable.myArray = [];  
myVariable.myArray[0] = "a";  
myVariable.myArray[1] = "b";
```



```
myVariable =  
{  
  "myNumber": 5,  
  myString: "bla",  
  myArray: ["a", "b", "c"]  
};
```

3 built-in plot-types

Line


Matrix

Polar

Adaptive procedure



Experiment Results



Parameter values: 54, 59, 64, 69, 71, 73, 75, 76, 77, 76, 77, 78, 77, 78, 79, 78, 77, 78, 79, 78, 79, 80, 79, 78, 79

Last value: 79


Reversals: 77, 76, 78, 77, 79, 77, 79, 78, 80, 78

**Mean (std)
last 6 reversals:** 78.5 (±0.957)

**Mean (std)
last 6 trials:** 78.8333 (±0.687)

Constant Procedure

Experiment Results




Correct % per stimulus									
B	Ch	D	F	M	N	P	S	T	V
25.00%	100.00%	100.00%	12.50%	75.00%	62.50%	100.00%	25.00%	25.00%	37.50%
Total percentage correct: 56.2500%									

	b	ch	d	f	m	n	p	s	t	v	w	z
b	2	0	0	2	0	0	0	0	0	1	0	0


ch	0	8	0	0	0	0	0	0	0	0	0	0
d	2	0	8	1	0	0	0	1	6	0	0	2
f	0	0	0	1	0	0	0	0	0	1	0	0
m	0	0	0	0	6	3	0	0	0	1	0	0
n	0	0	0	0	2	5	0	0	0	1	0	0
p	2	0	0	3	0	0	8	0	0	0	0	1
s	0	0	0	1	0	0	0	2	0	0	0	3
t	0	0	0	0	0	0	0	4	2	0	0	1
v	1	0	0	0	0	0	0	0	0	3	0	0
w	1	0	0	0	0	0	0	0	0	1	8	0
z	0	0	0	0	0	0	0	1	0	0	0	1

Trial	Stimulus	Answer	Correct
1	z	t	Incorrect
2	b	d	Incorrect
3	p	p	Correct
4	f	d	Incorrect
5	z	s	Incorrect
6	m	m	Correct
7	w	w	Correct



Localization


% correct per direction



Correct % per stimulus											
-90	-75	-60	-45	-30	-15	0	15	30	45	60	75
33.33%	0.00%	0.00%	33.33%	66.67%	0.00%	33.33%	33.33%	33.33%	0.00%	0.00%	33.33%
Total percentage correct: 25.6410%											



Overview of Resultsviewer internals



Overview of Resultsviewer internals

```
plots = {};
```

```
plots.types = ["line", "matrix", "polar", "text"];
```

```
plots.line = { data: [],  
               highChart: true,  
               prepare: prepareLine,  
               setData: null  
             };
```

Overview of Resultsviewer internals

```
plots = {};
```

```
plots.types = ["line", "matrix", "polar", "text"];
```

```
plots.line = { data: [],
```

```
  highChart: true,
```

```
  prepare: prepareLine,
```

```
  setData: null
```

```
};
```

```
plot()
```



```
preparePlotData()
```



```
setData()
```



Prepared plot data

prepare function generates data object

Line

```
plots.line.data =  
[ {  
  values: [],  
  reversals: [],  
  meanrevs: number,  
  meanrevstd: number,  
  meantrials: number,  
  meantrialstd: number  
}]
```

Matrix

```
plots.matrix.data =  
[ {  
  values: [ { x: xlabel ref,  
              y: ylabel ref,  
              z: x&y pairings  
            } ],  
  percentages: n,  
  xlabels: [],  
  ylabels: [],  
  raw: [][]  
}]
```

Polar

```
plots.polar.data =  
[ {  
  values: []  
}]
```

Setdata

Is only required if plot needs to be reconfigured/redrawn based on data

Matrix

```
function setDataMatrix(target, data, index)
{
    // x & y axis labels update
    target.xAxis[0].update({
        categories: data.xlabels,
        ceiling: data.xlabels.length - 1
    });
    target.yAxis[0].update({
        categories: data.ylabels,
        max: data.ylabels.length - 1,
        ceiling: data.ylabels.length - 1
    });
    return true;
}
```

Customizing with config object

Global

```
global =  
{  
  answers: ["pot", "pan", "kan", "kom"],  
  stimuli: [],  
  removefromanswer: ["word"],  
  multiprocedure: true/false  
}
```

Line

```
line =  
{  
  show: true/false,  
  parametername: "gain",  
  parameterunit: "dB",  
  parameterkey: "parametervalue",  
  trialsformean: n,  
  reversalsformean: n  
},
```

Customizing with config object

Matrix

```
matrix:  
{  
  show: true/false,  
  removefrommatrix: ["up", "down"]  
},
```

Polar

```
polar:  
{  
  show: true/false,  
  mindegree: n,  
  maxdegree: n  
},
```

How to customize

Inside .apx file

```
<results>  
  <page>resultsviewer.html</page>  
  <resultparameters>  
    <parameter name="line_reversalsformean">6</parameter>  
  </resultparameters>  
</results>
```



Format:

`name="plottype_parametername"`

for arrays such as global_stimuli

use , or | to input multiple values

e.g. pot,pan,kan,kom

How to customize

Custom resultsviewer

```
<results>  
  <page>my-resultsviewer.html</page>  
</results>
```

my-resultsviewer.html

```
...  
<script>  
  config.global.stimuli = ["-90","-45","0","45","90"]; //force list of possible stimuli  
  config.polar.show = true; //force polar plot to show  
</script>  
...
```

How to customize

In global config javascript file

Apex/resultsviewer/resultsviewer-config.js

- > Shows all possible configuration options
- > Will be applied to EVERY experiment/resultsviewer
- > Not recommended to edit this unless it's a specialized Apex install for one experiment

How to customize with js hooks

These are all shown in resultsviewer-config.js

> define these in your custom resultsviewer.html <script> tags

`mapAnswer(xmlstring)`

For custom XML answer-structures (when using pluginprocedures and such)

`resultsFilter(t)`

Filters out Trial objects based on certain criteria (function returns true to keep, false to remove)

`lineDataFilter(t)`

Same but specifically for Line plots

`matrixDataFilter(t)`

Same

`convertParameterValue(v)`

Automatically convert parametervalue with a formula

How to customize with special buttons

These are all shown in resultsviewer-config.js

> define these in your custom resultsviewer.html <script> tags

```
customButtons.push({
  name: "",
  label: "",
  replot: true/false,
  behavior: function()
  {
    //code
  }
});
```

Example

```
customButtons.push({
  name: "btn1",
  label: "Add series",
  replot: true,
  behavior: function()
  {
    var target = getChart("line");
    if (typeof target === "null")
      alert("Enable line first");

    target.addSeries( { data: [1,2,3,4,5,6] } );
    target.redraw();
  }
});
```

Advanced: Adding a new plottype

Highcharts.com

pick a plot, give it a name and add it to the plots object, eg "scatter"

```
plots.types.push("scatter");  
plots.scatter =  
{  
  data: [],  
  highChart: true,  
  prepare: function,  
  setData: function,  
  chartConfig: {}  
};
```